

Ior® Prestained Protein Marker

25 μL / 500 μL

Reference OZYC005-25 / OZYC005-500

SIZE: $25 \mu L / 500 \mu L$

STORAGE: Stable for up to 2 weeks at 25 °C or up to 3 months at 4°C. Store at -20°C for

long term storage.

SHELF-LIFE: 30 months

Product Description

The **ECOlor® Prestained Protein Marker** includes 10 prestained proteins of 3 colors. It covers a large range of molecular weights from 6.5 kDa to 270 kDa. When running on Tris-Glycine buffered SDS-PAGE, seven out of ten proteins in **ECOlor**® display **blue** bands since they are covalently labeled with a blue chromophore, the other three proteins show two orange bands (30 kDa and 270 KDa) and one green band (52 kDa). The ECOlor® Prestained Protein **Marker** is ideal for protein separation monitored on SDS-PAGE as it allows to verify the efficiency of Western transfer to any membrane (e.g. PVDF, nitrocellulose or nylon) as well as to estimate the proteins sizes. The **ECOlor® Pretained Protein Marker** is ready-to-use and includes the gel loading buffer.

Instruction for use

- 3 μL or 5 μL per loading for clear visualization during electrophoroesis on 15-well or 10-well mini-gel, respectively
- ✓ 2.5 µL per well for general Western transferring.

Storage Buffer

62.5 mM Tris-H₃PO₄ (pH 7.5 at 25 °C), 1 mM EDTA, 2% SDS, 10 mM DTT, 1 mM NaN₃, 33% glycerol.

Quality Control


Under suggested conditions, **ECOlor® Prestained Protein Marker** resolves 10 major bands in 4 - 12% Bis-Tris Gel (MES buffer) and after Western blotting to the nitrocellulose membrane.

Technical support Email: tech@ozyme.fr Phone: +33 (0) 1 34 60 60 24 Order information Email: commande@ozyme.fr

Phone: +33 (0) 1 34 60 15 16 Fax: +33 (0) 1 34 60 92 12

Guide for Molecular Weight Estimation* (kDa)

Migration patterns of ECOlor® Prestained Protein Marker in different electrophoresis conditions are showed below:



*5 µL of ECOlor® Prestained Protein Marker per lane

Note: The apparent molecular weight (kDa) of each protein has been determined by calibration against unstained protein standards; supplemental data should be considered for more accurate adjustment in different electrophoresis conditions.

Ordering Information

OZYME

6 bd Georges Marie Guynemer

ZAC Charles Renard - Bât. G

78210 SAINT-CYR-L'ÉCOLE - France

Email: commande@ozyme.fr

Phone: +33(0) 1 34 60 15 16 Fax: +33 (0) 1 34 60 92 12

Document version

OZYC005-0120

PRODUCT FOR RESEARCH USE ONLY

Techhical support Email: tech@ozyme.fr Phone : +33 (0) 1 34 60 60 24 Order information Email: commande@ozyme.fr Phone: +33 (0) 1 34 60 15 16 Fax: +33 (0) 1 34 60 92 12