


# Azure Biosystems IQOQ

## Validation package for Cielo Real-Time PCR Systems

The Azure Biosystems® IQOQ package is an essential tool for verifying any Cielo™ Real-Time PCR System. The validation package contains the documentation required for regulatory compliance, SOPs, and laboratory management. It includes the protocols, checklists, and certificates necessary to satisfy your validation requirements.


### FEATURES

- Detailed procedures for Installation Qualification (IQ) and Operational Qualification (OQ)
- For use with any Cielo Real-Time PCR system
- Easy-to-use checklists with clear instructions



### IQOQ Protocol

Includes the procedures for Installation Qualification and Operational Qualification



### IQOQ Checklist & Logbook

Includes checklists, and document control logs



### IQOQ Certificates



### OQ Validation Plate Reagents

(sold separately)

# IQOQ phases for Cielo Real-Time PCR Systems Models 3 and 6

## 1. Installation Qualification (IQ)

During IQ phase, the equipment is inspected to ensure it meets the requirements of Azure Biosystems and has been installed properly.

IQ Validation includes:


- Examining equipment design and supplied documentation
- Assuring that the work area has the required space for the system to operate at peak performance
- Determining installation parameters
- Identifying required environmental controls and procedures
- Installing the device
- Verifying correct installment of the device
- Establishing the procedures for monitoring, operating, and control of the device including responsible operators

## 2. Operational Qualification (OQ)

OQ phase tests the newly installed device to ensure it operates as expected and gauges the conditions under which it will continue to do so. Test runs will determine that the system is operating as expected. If OQ is successful, it will validate that the device is operating at optimal performance.

OQ Validation includes:

- Validating system internal diagnostics
- Validating application operations



## Ordering Information

Part Number	Name	Description
AIQ317	IQOQ	IQOQ procedures for the Cielo models 3 and 6. Includes procedural manual, checklist, and certifications.
AIQ316	IQOQ Service	IQOQ Service performed by Azure Biosystems Personnel
AIQ318	OQ Validation Plate Reagents	Reagents required for validation plate for OQ procedure

**Nous contacter**