

OZYME

Des femmes et des hommes
au service de vos recherches

COMPATIBILITY CHART OF PIPETTE TIPS

We confirm that pipette tips fit on following pipette

Brand Pipettes Fix/Adjustable	Pipette tips with different liquid volume				
Eppendorf	10µl	20µl	100µl	200µl	1000µl
10µl	√				
20µl		√			
100µl		√	√		
200µl		√	√	√	
1000µl					√
Dragon	10µl	20µl	100µl	200µl	1000µl
10µl	√				
20µl		√			
100µl		√	√		
200µl		√	√	√	
1000µl					√
Gilson	10µl	20µl	100µl	200µl	1000µl
10µl	√				
20µl		√			
100µl		√	√		
200µl		√	√	√	
1000µl					√
Rainin	10µl	20µl	100µl	200µl	1000µl
10µl	√				
20µl		√			
100µl		√	√		
200µl		√	√	√	
1000µl					√
Thermo Finnpipette	10µl	20µl	100µl	200µl	1000µl
10µl	√				
20µl		√			
100µl		√	√		
200µl		√	√	√	
1000µl					√
Remarks:					
1. These pipette tips are compatible with brands pipettes including <i>Eppendorf, Dragon, Gilson, Rainin, Thermo, Socorex, Transferpette</i> for single and multi-channel pipettes, except for <i>Rainin and Socorex</i> (single-channel pipettes only)					
2. For other brands pipettes incl <i>Sartorius, Biohit</i> , we will list them based on more testing.					

Nous contacter

Service client - commande : commande@ozyme.fr

Service technique : Réactifs : 01 34 60 60 24 - tech@ozyme.fr - Instrumentation : 01 30 85 92 88 - instrum@ozyme.fr