

USAGE LIMITS - BEAD BEATING TUBES OZYME

Reference	Name	Tube Volume	Dry Grinding		Liquid Grinding		
			Max Speed	Max Cycle	Max Speed	Max Cycle	Max Volume
OZYE001-100	Ozyme Lysing Kit B, 0,5 mm glass beads (micro-organism lysing)	2 mL	10 000 rpm	10 x 30 s	10 000 rpm	10 x 30 s	1,2 mL
OZYE002-100	Ozyme Lysing Kit C, 2,8 mm ceramic beads (hard tissue homogenizing)		8 800 rpm	10 x 30 s	10 000 rpm	10 x 30 s	1,2 mL
OZYE003-100	Ozyme Lysing Kit D, 1,4 mm ceramic beads (soft tissue homogenizing)		10 000 rpm	10 x 30 s	10 000 rpm	10 x 30 s	1,2 mL
OZYE004-100	Ozyme Lysing Kit M (Mix), 1,4 mm & 2,8 mm ceramic beads (tissue homogenizing)		8 800 rpm	10 x 30 s	10 000 rpm	10 x 30 s	1,2 mL

OZYME - 6 Bd Georges Marie Guynemer ZAC CHARLES RENARD - Bât.G - 78210 SAINT-CYR-L'ECOLE

Tél. : 01 34 60 24 24 - Fax : 01 34 60 92 12 - www.ozyme.fr/info

Service client - commande : commande@ozyme.fr

Service technique : Réactifs : 01 34 60 60 24 - tech@ozyme.fr / **Instrumentation** : 01 30 85 92 88 - instrum@ozyme.fr