

Protein Safe Lysing kits



Available for CK-14 & CK28-R lysi

Protein Safe Lysing kits

Prevent protein degrading by shutting down unwanted protease activity & ensure your protein stays intact for your experiment

Protein Safe is efficient against :

- Serine Proteases
- Aminopeptidase B & Leucine Aminopeptidase
- Cysteine Proteases & Trypsin-like Proteases
- Aspartic Proteases

Product references

- Protein Safe Hard tissue homogenizing CK28-R: P000972-LYSKO-A.0
- Protein Safe Soft tissue homogenizing CK14: P000973-LYSKO-A.0

PRECELLYS LYSING KITS RANGE




- 5** types of tube volumes
(0,3mL in well-plate format, 0,5, 2, 7 & 15mL)
- 4** grinding materials
(ceramic, glass, garnet and stainless steel)
- 30+** lysing kits

COMPOSITION:

AEBSF, aprotinin, bestatin, E-64, Leupeptin, Pepstatin A

Discover our Life Sciences product range



1500+ SCIENTIFIC DOCUMENTS

SIGN UP TO THE APPLICATION CENTER
www.bertin-technologies.com



NOUS CONTACTER

VOIR LES ÉQUIPEMENTS

OZYME
Des femmes et des hommes
au service de vos recherches

Service client - commande: commande@ozyme.fr

Service technique:

Réactifs: [01 34 60 60 24 - tech@ozyme.fr](mailto:01_34_60_60_24_tech@ozyme.fr)

Instrumentation: [01 30 85 92 88 - instrum@ozyme.fr](mailto:01_30_85_92_88_instrum@ozyme.fr)